# FRANKA PRODUCTION 3 FRANKA RESEARCH 3

---

## Quick Guide for Installation

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## NOTICE

This Guide is not the manual of your Franka robot and does not replace it.

## INDICATIONS

---

Before operating robots from Franka Robotics, all persons working with them have to read the corresponding manual carefully. Read or download the manual (translations available) here: [www.franka.de/documents](http://www.franka.de/documents).

Don't hesitate to contact your sales partner or integrator if you need more detailed information. Alternatively, you can get in touch with us directly by sending an e-mail to: [support@franka.de](mailto:support@franka.de).

## NOTICE

Please carefully follow all instruction, warnings and other notes in this document and the official manual. Any deviation from the steps prescribed in these documents constitutes a misuse (as described in the paragraph "Misuse" of the corresponding manual). Franka Robotics does not assume any warranty and/or liability for any possible consequences.

Get your manual and additional supporting material  
in English and other languages.  
[www.franka.de/documents](http://www.franka.de/documents)

## QUICK GUIDE FOR INSTALLATION

**This booklet will guide you  
through the product manual's steps  
to get your system up and running.**

Document number: R02040  
Release Version: 1.0

The English documentation is the ORIGINAL DOCUMENTATION.  
Other languages are translations of the original document.

See relevant product manual at paragraph:  
Equipment Overview


**⚠ DANGER!**

Safety-related aspects must be considered when planning and designing the application and performing the hazard & risk assessment for completed machinery.


## System components and connections

1	Arm
2	Emergency Stop Device
3	Connecting cable
4	Control
5	Ethernet (Network)
6	Power switch
7	Power cable
8	Mains power outlet
9	Interface device ( <i>not included</i> ) with Franka UI
10	External Enabling Device
11	Emergency Unlocking System
12	Connection to functional earth


# Overview


Left view

Front view

Back view

See relevant product manual at paragraph:  
Unpacking the Equipment

### WARNING

#### Heavy equipment

Due to the dead weight and partly due to the geometric design, lifting and handling the equipment may cause back injury and, if it falls, serious injury to fingers, hands, toes, and feet.

- When transporting the equipment, always wear personal protective equipment (e.g., safety shoes).
- Always lift the equipment with the help of a second person.
- Follow the existing company regulations on lifting loads and personal protective equipment.

### NOTICE

Always keep the original packaging in case of relocating the robot.

### NOTICE

Franka Hand is not part of the certified machinery.


## Unpacking

#### Procedure:

- Take off the top lid of the transport box (A).
- Lift the top inner boxes (i.e., 2-5) and put them aside (B).
- Pull the Arm box (i.e., 1) out of the lower transport box (C).

1	Arm
2	Control
3	Connecting cable
4	Emergency Stop Device & External Enabling Device
5	Optional (e.g., Franka Hand)

# Unpacking

See relevant product manual at paragraph:  
Preparing the Installation Site

#### NOTICE

The heat produced by power electronic components and modules inside the Arm is dissipated via the Arm's surface. The heat produced by power electronic components and modules inside the Control is dissipated through an internal ventilation system.

- Install the Arm and the Control at a place in an adequately ventilated place.
- Do not expose the Arm and the Control to direct sunlight.
- Do not repaint, paste up, wrap up the Arm.
- Place the Control at a sufficient distance between the front/back fans and covering components (40 mm on both sides).
- Make sure that the Control's fans are not covered with dirt.

#### NOTICE

Install the Arm in a position that allows for ergonomic human-robot interaction.

## Preparing the installation site

#### Required material:


- Detailed technical drawing for baseplate.

#### Procedure:


- Use the technical drawing for positioning and creating the threads and fitting holes.

The Arm is equipped with highly sensitive sensor technology and fine-tuned control algorithms, which require installation on a stable, leveled ( $\leq 0.1^\circ$  tilt angle), non-moving, and non-vibrating platform in an upright position. Max. forces (i.e., 5-8) must be supported during static and dynamic operation.

# Step 1

Print 1:1 to use as drilling template  
[www.franka.de/documents](http://www.franka.de/documents)

1	Front	
2	Holes for screws (M8)	
3	Thread for functional earth (M5)	
4	Holes for alignment pins (Ø6H7)	
5	Tilting torque	280 Nm
6	Torque around axis	190 Nm
7	Horizontal force	300 N
8	Vertical force	410 N
9	Leveled surface	≤ 0.1° tilt angle

See relevant product manual at paragraphs:  
Unpacking the Equipment, Positioning the Control

### WARNING!

#### Heavy equipment

Due to the dead weight and partly due to the geometric design, lifting and handling the equipment may cause back injury and, if it falls, serious injury to fingers, hands, toes, and feet.

- When transporting, mounting, or demounting the equipment, always wear personal protective equipment (e.g., safety shoes).
- Always lift the equipment with the help of a second person.
- The equipment must be placed on even surfaces to prevent it from tilting or sliding.
- Follow the existing company regulations on lifting loads and personal protective equipment.




### NOTICE

Allowed mains frequency: 50 - 60 Hz  
Supply voltage: 100 - 240 VAC

## Unpacking and placing the Control

- Carefully open the box by removing the sealing adhesive strips on top of the cardboard box.
- Open the foil coating.
- Remove the power cable and top lid (A).
- Remove the top protective layer (B).
- Grasp the Control at the indicated grasping positions and carefully lift it out of the bottom protective layer and set it aside (C).
- Remove the foam packaging from the Control.
- Place the Control horizontally in its designated position, or attach it to a rack designed for 19" units.

## Step 2

See relevant product manual at paragraph:  
**Unpacking the Equipment**

### **WARNING!**

#### **Heavy equipment**



Due to the dead weight and partly due to the geometric design, lifting and handling the equipment may cause back injury and, if it falls, serious injury to fingers, hands, toes, and feet.

- When transporting, mounting, or demounting the equipment, always wear personal protective equipment (e.g., safety shoes).
- Always lift the equipment with the help of a second person.
- The equipment must be placed on even surfaces to prevent it from tilting or sliding.
- Follow the existing company regulations on lifting loads and personal protective equipment.

## **Unpacking and lifting the Arm**

Always lift the Arm in the positions intended for lifting, as not to overstress the the Arm joints during handling and lifting. In particular, the Arm may never be carried in the extended position with one person holding each end of the Arm.

### Step 3

See relevant product manual at paragraphs:  
**Mounting the Arm, Wiring and Electrical Installation**

## NOTICE

Ensure that the maximum forces and torques are supported by the mounting surface during static and dynamic operation. For more information, see Step 1.

## Mounting the Arm

The Arm must be connected securely to the baseplate with 4 suitable screws. For this purpose, 4 drill holes with a diameter of 9 mm are provided in the base flange of the Arm.

### Required material:

- Washers and screws depend on the surface the robot is mounted on. Please check the table in paragraph “Mounting the Arm” for details.
- 1x cylindrical head screw with hexagon socket M5x8 (strength class 8.8 A2K)
- 1x tooth washer M5 (strength class A2K)
- Torque spanners to tighten screws with 30 Nm


### Precondition:

- Two people for placing the Arm
- Prepared baseplate

### Procedure:

- Lift the Arm.
- Carry the Arm to its designated position.
- Align the Arm accordingly to the pre-fixed holes on the baseplate.
- Person 1: Hold the Arm.  
Person 2: Use the four screws to mount it to the baseplate.
- Connect the functional earth to the base of the Arm. We recommend using a minimum 1.5 mm<sup>2</sup> section copper (Cu) cable with a maximum length of 5 m.
- Make sure that the Emergency Unlocking Tool is inserted into its holder. For more information on its use, see paragraph “Manually Moving the Arm”.

## Step 4

See relevant product manual at paragraph:  
**Wiring and Electrical Installation**



## NOTICE

Only Connecting cables provided by Franka Robotics shall be used for electrical connection between Arm and Control.

## Plugging in the Connecting cable

- Carefully place the cable's plug onto the connector X1, with the triangular marking pointing upward. Make sure that you use the correct pins/holes coupling (A).
- The plug itself is pulled into the connector by turning the movable front part of the plug. Turn hand-tight (B).
- Apply the same principle to connect the other end of the Connecting cable with the C1 connector on the front of the Control (C, D).

## Step 5

See relevant product manual at paragraph:  
**Wiring and Electrical Installation**

## INDICATIONS

---

Emergency stop devices need to be installed according to the generally valid and accepted engineering standards, e.g., European standards EN 60204 and related. The Emergency Stop Device provided by Franka Robotics needs to be connected to the X3 port. Other devices than the Emergency Stop Device provided by Franka Robotics may be connected to the X3 port as well.



The devices connected to the emergency stop signal must follow EN 60947-5-5 or EN 62061. Connected emergency stop devices must be placed such that they are easily reachable in case of emergency. Disconnected emergency devices must be removed to prevent misunderstandings.

For more information on wiring, see “Wiring for X3” in paragraph “Safety Settings and Watchman”.

## Completing the installation

- Connect the External Enabling Device to the connector X4 on the base of the Arm.
- Ensure that the guide pins are pointing in the right direction (A).
- The plug itself is pulled into the connector by turning the movable front part of the plug. Turn hand-tight (B).
- Connect the Emergency Stop Device to the connector X3 on the base of the Arm.
- Ensure that the guide pins are pointing in the right direction (C).
- By turning the movable front part of the plug, the plug itself is pulled into the connector. Turn hand-tight (D).
- Place both devices within easy reach of the user.

## Step 6

See relevant product manual at paragraphs:  
[Wiring and Electrical Inst., Switching On, Correct Inst. Site](#)

**⚠ DANGER!**

**Electrical shock**

Short circuit due to developed condensation when the device is transported from colder to warmer and more humid surroundings. Risk of life-threatening injuries due to electrical shock.

- Leave devices to acclimatize after transport.
- Do not switch on wet devices.

**⚠ DANGER!**

**Damaged wires or inadequate electrical installation**

Risk of personal injury by electric shock, as well as material damage.

- Only use the system in sound technical condition.
- Only have the emergency stop system installed by qualified personnel.
- Check cables and electrical installations.

**⚠ CAUTION!**

**Exposed wires and cables**

Operators may trip and fall due to exposed wires and cables in the working space.

- Always lay cables safely.

## Switching on the system




**Precondition:**

- Cables need to be plugged in correctly.
- The external power supply must be connected (A).
- The knob of the Emergency Stop Device needs to be pulled up (B).
- Leave the maximum space (see paragraph “Correct Installation Site” in product manual).

**Procedure:**

- Switch on the Control. The status lights will flash white (C).

## Step 7

See relevant product manual at paragraphs:  
**Connecting a User Interface Device, First Start**

### **DANGER!**

Safety-related aspects must be considered when planning and designing the application and performing the hazard & risk assessment for completed machinery.

### **NOTICE**

Do not place devices (e.g., laptop) in the hazardous zone of the robot.

### **INDICATIONS**

Register to Franka World – <https://franka.world/> – to keep your system up to date with the latest software, explore and purchase a steadily growing portfolio of products, and access documentation, tutorials, code, and tools.

## **Configuring the system**

### **Precondition:**

- Access to a browser, e.g., Chrome, Chromium, or Firefox.

### **Procedure:**



- To start the initial configuration, connect the interface device to the X5 connector on the Arm base via Ethernet cable (*not included*).
- Connect (optional) the Control to the network via ethernet cable (*not included*).
- Go to your preferred web browser.
- Enter the following URL: **robot.franka.de**
- The First Start page will appear.
- From here, the product manual can be accessed, even offline.

**Before continuing to work with the robot, read the product manual carefully.**

- Go through the steps of the First Start configuration, as described in the chapter “Operation” of the product manual.
- After rebooting the system and creating an administrator role, Desk is displayed in the web browser. The status lights will glow solid blue.
- Relevant safety settings are necessary for unrestricted operation. Follow required steps, as described in chapter “Watchman” of the product manual.

## Final Step

Enjoy using your Franka robot! We look forward to getting to know about your success! Share your experience with us, and we will amplify the reach on our own digital and social media platforms.

Just tag **@frankarobotics** or **#frankarobotics** in your post, and we will ensure your endeavors will get international visibility!

For more marketing and PR opportunities, reach out to **marketing@franka.de**.

Print for use as a quick reference  
[www.franka.de/documents](http://www.franka.de/documents)

## Color code of status lights

Color	Category
White	Not active; ready to start execution
White (slow flashing)	Not active; booting or shutting down
White (fast flashing)	Not active; updating
Blue	Brakes engaged
Blue	Ready to execute task
Blue (slow flashing)	Opening brakes
Blue (slow flashing)	Collaborative operation (no task active)
Green	Automatic execution
Green (slow flashing)	Collaborative operation (task active)
Green (fast flashing)	Attention: Automatic execution starting
Yellow	Warning
Yellow (slow flashing)	Warning; user interaction needed
Pink (slow flashing)	Conflicting inputs
Red	Error
Red (slow flashing)	Safety violation / application error

Slow flashing: 0.6 Hz, twice per 3 seconds  
Fast flashing: 2 Hz, twice per second

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**Franka Robotics GmbH**